

# **Affinity-based Dispersive Pipette Extraction for Automated Purification**

**Todd Mullis<sup>1</sup>**, Sunil Hwang<sup>1</sup>, Michael Mouradian<sup>2</sup>, L. Andrew Lee<sup>\*1</sup> <sup>1</sup>IMCS, LLC - Irmo, SC; <sup>2</sup>Hamilton Robotics - Reno, NV





- 20 µL of 50% MabSelect SuRe<sup>TM</sup> LX slurry
- Binding buffer: 20 mM sodium phosphate, 0.15 M sodium chloride
- Elution buffer: 0.1 M sodium citrate, pH 3.5
- Neutralization buffer: 1 M tris, pH 9.0

### IMCS cobalt-IMAC tips on automatic pipetting system

- 1 mg of cell lysate
- 50 µL of 50% cobalt-IMAC slurry
- Binding buffer: 20 mM sodium phosphate, 0.5 M NaCl, 20 to 40 mM imidazole, pH 7.4
- Elution buffer: 20 mM sodium phosphate, 0.5 M NaCl, 500 mM imidazole, pH 7.4

### RESULTS

Immunoglobulins from human serum samples were purified using IMCStips loaded with resin modified with Protein A on multi-channel liquid handling system (Nimbus and VIAFLO 96) within 30 minutes from start to finish. The extraction process in comparison to conventional spin column format, did not require additional incubation time and recoveries were comparable or better than conventional formats (*Figure 3*). His-tagged enzyme was also purified from crude cell lysates (Figure 6) in similar fashion on both multi-channel liquid handling systems within 30 minutes using IMAC resin. The final purity of the target enzymes (> 80%) are shown in Figure 4b. Further optimization of target proteins could be achieved by modifying the wash buffer compositions to achieve higher purity.



The workflow demonstrated here show the potential to achieve 96 sample purification within 30 minutes by leveraging dispersive pipette extraction on multi-channel liquid handling systems.






Figure 3. Side by side comparison of Protein A IMCStip elution and Protein A spin column. 20  $\mu L$  of resin slurry was used for  $10 \mu L$  of human serum




Figure 5. A comparison of protein recovery for 24 serum samples enriched using the 300 µL Protein A IMCStip. These samples were biological replicates to show the reproducibility of the extraction. The calculated CV value for this sample set was 6.5%




Figure 6. Before and after enrichment of 8 different recombinant his-tagged proteins in lysate using the cobalt-IMAC IMCStips.

## CONCLUSIONS

We successfully developed a fully automated high-throughput affinity based protein extraction from complex matrices. These findings will allow for 30 minute complete sample preparation leading faster screening and downstream applications.

Next steps are to use different types of affinity resins in order to purify different targets in complex matrixes.

## REFERENCES

- 1. GE Healthcare. (2015). Antibody Purification. GE Healthcare Bio-Sciences AB: Uppsala, Sweden. https://proteins.gelifesciences.com/~/media/protein-purificationib/documents/handbooks/antibody\_purification\_handbook.pdf?la=en
- 2. GE Healthcare. (2016). Affinity Chromatography: Vol. 2: Tagged Proteins. GE Healthcare Bio-Sciences AB: Uppsala, Sweden. http://proteins.gelifesciences.com/~/media/protein-purification-ib/documents/handbooks/affinity-chromatography-handbook-vol-2.pdf
- <sup>1</sup>Contact: Andrew Lee lee@imcstips.com
- <sup>2</sup>Contact: Michael Mouradian michael.mouradian@hamiltoncompany.com
- Artwork and layout created by Carmen Adamson.

© 2017 IMCS, LLC. All rights reserved. IMCStips<sup>™</sup> is a trademark of Integrated Micro-Chromatography Systems, LLC. Microlab<sup>®</sup> NIMBUS <sup>®</sup> is a registered trademark of Hamilton Robotics (Reno, NV)



